



Email : sales@laboratoryinstrumentindia.com

Product Name :
Flexural And Transverse Multipurpose Testing Machine

Product Code :
LBNY-0005-1220090



Description :

Load is measured by a high accuracy electric strain cell, eliminating the pistons weight and friction.
Multipurpose high flexibility testing machine.
C-shaped open frame for an easy and fast positioning of the specimen between the rollers.
The frame is closed by a hydraulic vertical rod, granting high rigidity.
Maximum vertical daylight between upper/lower rollers: 263 mm
Horizontal clearance (between uprights): 1040 mm
Maximum capacity: 320 KN
Ram travel: 110 mm
Graduated scales are foreseen to get easy rollers adjustment.
Simple action piston with counterweights to optimize frictions.
Possibility to easily place in the centre one upper roller for centre point loading.

Technical Specification :

Power supply: 230V 1ph 50Hz 750W
Calibration accuracy: grade 1.0
Frame dimensions: 1700 x 1470 x 1557 mm
Frame weight: 800 kg.

Laboratory Instrument India

