



Email : sales@laboratoryinstrumentindia.com

Product Name :
Flexural Testing Machine 200 KN Digital Model High
Stability

Product Code :
LBNY-0005-1220091



Description :

This two frames Group offers the considerable advantage to perform compression tests on concrete cube, cylinder and block specimens flexural tests on concrete beams and by using suitable accessories, to perform also the following tests: Compression on portions of 40.1x40x160 mm mortar broken in flexure. Flexural high stiffness frame 200 KN capacity complete with pressure transducer, used in conjunction with a digital compression machine.

Technical Specification :

Compression on 70 mm mortar cubes,
Compression on 50 mm mortar cubes
Splitting tensile on concrete cubes and concrete block pavers
Splitting tensile on cylindrical specimens dia. 100, 150, 160 mm.



Laboratory Instrument India