



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Concrete Compression Machine 2000 KN High Stability

**Product Code :**  
LBNY-0005-12200094



**Description :**

Compression testing machine 2000 KN motorized with control unit high stability to test blocks max. 500x300 mm, cubes up to 200 mm side and cylinders up to dia. 160x320 mm.

**Technical Specification :**

Concrete Compression Machine 2000 KN High Stability



**Laboratory Instrument India**