



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Concrete Compression Machine 2000 KN Motorized

**Product Code :**  
LBNY-0005-1220095



**Description :**

Compression testing machine 2000 KN capacity motorized with Servo-Plus evolution control unit to test blocks max. 500x300 mm, cubes up to 300 mm side and cylinders up to dia. 160x320 mm.

**Technical Specification :**

Concrete Compression Machine 2000 KN Motorized



**Laboratory Instrument India**