



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Concrete Compression Machine 4000 KN Motorized High  
Stability Semi Automatic Model

**Product Code :**  
LBNY-0005-12200096



**Description :**

Compression testing machines motorized with control unit (semi-automatic model) to test cubes up to 200 mm side and cylinders up to dia. 250x500 mm.

**Technical Specification :**

Concrete Compression Machine 4000 KN Motorized High Stability Semi Automatic Model



**Laboratory Instrument India**