



Email : sales@laboratoryinstrumentindia.com

Product Name :
Concrete Compression Machines 3000 KN Motorized To
Test Cubes Cylinders And Blocks

Product Code :
LBNY-0005-1220097



Description :

Compression machine 3000 KN motorized with control unit, to test blocks max. 500x300 mm, cubes up to 300 mm side and cylinders up to dia. 160x320 mm

Technical Specification :

Concrete Compression Machines 3000 KN Motorized To Test Cubes Cylinders And Blocks



Laboratory Instrument India