



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Concrete Compression Machines 3000 KN 1 Gauge To  
Test Cubes And Cylinders

**Product Code :**  
LBNY-0005-1220098



**Description :**

Compression machine 3000 KN motorized 1 gauge to test cubes up to 200 mm side and cylinders up to dia. 160x320 mm.

**Technical Specification :**

Concrete Compression Machines 3000 KN 1 Gauge To Test Cubes And Cylinders



**Laboratory Instrument India**