



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Concrete Compression Machine 2000 KN Motorized Semi Automatic Model

**Product Code :**  
LBNY-0005-122000100



**Description :**

Compression testing machine 2000 kN capacity motorized with control unit (semi-automatic model), to test cylinders up to dia. 160x320 mm and cubes up to 150 mm side.

**Technical Specification :**

Concrete Compression Machine 2000 KN Motorized Semi Automatic Model



**Laboratory Instrument India**