



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Concrete Compression Machines 3000 KN Motorized To  
Test Cubes And Cylinders

**Product Code :**  
LBNY-0005-122000101



**Description :**

Compression machine 3000 kN, motorized with control unit, to test cubes up to 200 mm side and cylinders up to dia. 160x320 mm.

**Technical Specification :**

Concrete Compression Machines 3000 KN Motorized To Test Cubes And Cylinders



**Laboratory Instrument India**