



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Concrete Compression Machine 3000 KN Motorized High  
Stability Semi Automatic Model

**Product Code :**  
LBNY-0005-122000102



**Description :**

Compression testing machines 3000 KN motorized with control unit semi-automatic model, high stability, to test cubes up to 200 mm side and cylinders up to dia. 160x320 mm.

**Technical Specification :**

Concrete Compression Machine 3000 KN Motorized High Stability Semi Automatic Model



**Laboratory Instrument India**