



Email : sales@laboratoryinstrumentindia.com

Product Name :
Concrete Compression Machine 3000 KN Motorized High
Stability Semi Automatic Model

Product Code :
LBNY-0005-122000102



Description :

Compression testing machines 3000 KN motorized with control unit semi-automatic model, high stability, to test cubes up to 200 mm side and cylinders up to dia. 160x320 mm.

Technical Specification :

Concrete Compression Machine 3000 KN Motorized High Stability Semi Automatic Model



Laboratory Instrument India