



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Concrete Compression Machine 1500 KN Motorized 1 Gauge

**Product Code :**  
LBNY-0005-122000103



**Description :**

Compression machine 1500 KN motorized 1 gauge, to test cubes up to 150 mm side and cylinders up to dia. 160x320

**Technical Specification :**

Concrete Compression Machine 1500 KN Motorized 1 Gauge



**Laboratory Instrument India**