



Email : sales@laboratoryinstrumentindia.com

Product Name :
Compression Machines 2000 KN To Test Cylinders And
Cubes High Stability Motorized

Product Code :
LBNY-0005-122000105



Description :

Compression testing machines 2000 KN Capacity tested for high stability to test cubes up to 200 mm side and cylinders up to dia. 160x320 mm.

Technical Specification :

Hydraulic device to stop the pistons stroke at its max excursion, to avoid pumping the piston out of the cylinder.
Gauges dia. 250 mm with specific resistance scales for cubes 150 mm and cylinders dia. 150 and 160 mm
Gauges divisions: 2000 KN div. 5 KN - 600 KN div. 2 KN
Compression platens dia. 287x60 mm
Max. vertical daylight: 336 mm
Power supply: 230 V 1 ph 50 Hz 750 W
Calibration accuracy: Grade 1.0
Max. Ram travel 60 mm approx.
Dimensions: 690x400x1400 mm
Weight: 850÷920 kg.



Laboratory Instrument India

