



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Concrete Compression Machine 3000 KN High Stability  
Motorized

**Product Code :**  
LBNY-0005-122000106



**Description :**

Compression machine 3000 KN high stability, motorized, with control unit to test blocks max. 500x300 mm, cubes up to 200 mm side and cylinders up to dia. 160x320 mm.

**Technical Specification :**

Concrete Compression Machine 3000 KN High Stability Motorized



**Laboratory Instrument India**