



Email : sales@laboratoryinstrumentindia.com

Product Name :
Compression Machines 2000 KN To Test Cubes

Product Code :
LBNY-0005-122000107



www.laboratoryinstrumentindia.com

Description :

Compression testing machines 2000 KN capacity to test cubes up to 150 mm side and cylinders up to dia. 160x320 mm.

Technical Specification :

Hydraulic device to stop the pistons stroke at its max excursion to avoid pumping the piston out of the cylinder.
Gauges dia. 250 mm with specific resistance scales for cubes 150 mm and cylinders dia. 150 and 160 mm
Gauges divisions: 2000 KN div. 5 KN - 600KN div. 2 KN
Max. vertical daylight: 336 mm
Compression platens dia. 216 mm
Max. Ram travel 55 mm approx.
Calibration accuracy: Grade 1.0
Power supply: 230 V 1 ph 50 Hz 750 W
Dimensions: 690x400x1320 mm
Weight: 700 Kg.



Laboratory Instrument India

