



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Universal Digital Tester With Microprocessor For Load Cells

**Product Code :**  
LBNY-0005-122000112



**Description :**

The instrument foresees three memorized cycle verification program composed of ten measurement each. This user friendly menu driven digital display, connected to load cells allows to perform an accuracys verification of the loads measured from machines under control and it allows to produce the relative certificate.

At the end of the test the unit automatically elaborates the stored value and displays:

Average measured load

Accuracy in %

Effective applied load

Repeatability

Measured load.

**Technical Specification :**

Universal Digital Tester With Microprocessor For Load Cells



**Laboratory Instrument India**