



Email : sales@laboratoryinstrumentindia.com

Product Name :
Workshop Dye Penetrant

Product Code :
LBNY-0004-740030



Description :

Workshop Dye Penetrant

Technical Specification :

Dye penetrant
Standard for welding
Other Names: Dye Penetrant testing kit
Grade Standard: Industrial Grade
Application: Crack testing in welds, surface defect inspection on steel plates.

Workshop Dye Penetrant, Workshop Dye Penetrant Equipments, Workshop Dye Penetrant Tools, Workshop Dye Penetrant Tool Kits, Workshop Dye Penetrant Manufacturers, Workshop Dye Penetrant Suppliers from India, China, Kenya



Laboratory Instrument India