



Email : sales@laboratoryinstrumentindia.com

Product Name :
Poker Vibrator

Product Code :
LBNY-0005-122000114



Description :

Poker Vibrator, portable, petrol operated.
Motor, 4-stroke, 1,6HP, 35,8cc
Flexible shaft: 2 metres long.
Frequency: 10000 vibrations/min.
Tip dimensions: 25 mm dia. by 250 mm long.
Supplied complete with knapsack.

Technical Specification :

Weight: 7 kg



Laboratory Instrument India