



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Compression Machines For Cylinders And Cubes Touch  
Screen

**Product Code :**  
LBNY-0005-122000115



**Description :**

Compression testing machines 1500 kN capacity to test cubes up to 150 mm side and cylinders up to dia. 160x320 mm.

**Technical Specification :**

Hydraulic device to stop the pistons stroke at its max excursion to avoid pumping the piston out of the cylinder.

Max. vertical daylight: 336 mm

Max. ram travel 55 mm approx.

Power supply: 230 V 1 ph 50 Hz 750 W

Compression platens dia. 216 mm

Calibration accuracy: Grade 1.0

Dimensions: 630x350x1260 mm

Weight: 540÷580 kg.



**Laboratory Instrument India**

