



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Concrete Pipe Testing Machine Capacity Load

**Product Code :**  
LBNY-0005-122000117



**Description :**

Designed to test concrete sewer and drain pipes used in drainage works, water and irrigation supply systems etc.

The machine is composed of two parts:

Steel Testing Frame.

Electro-Hydraulic loading and control cabinet.

**Technical Specification :**

Concrete Pipe Testing Machine Capacity Load



**Laboratory Instrument India**