




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Splitting Tensile Test Device For Concrete Block Pavers

**Product Code :**  
LBNY-0005-122000120



**Description :**

Splitting Tensile Test Device For Concrete Block Pavers is directly fixed on the compression platens of the block testers having 2000 kN or 3000 kN capacity. Splitting tensile test device to perform tests on concrete block pavers having max. dimensions 300 x 500 mm, and for tests on concrete cube specimens 100, 150, 200 mm, and any type of block and prismatic specimens.

**Technical Specification :**

Splitting Tensile Test Device For Concrete Block Pavers



**Laboratory Instrument India**