



Email : sales@laboratoryinstrumentindia.com

Product Name :
Concrete Flexural And Transverse Machine 150 KN

Product Code :
LBNY-0005-122000123



Description :

Concrete Flexural And Transverse Machine 150 KN with control unit to perform flexural tests on concrete beam specimens max. size 200x200x800 mm, kerbs, tiles, slabs, masonry units, flat blocks, flagstones, pipes, and any type of material having max. size 550xh550 mm.

Technical Specification :

Concrete Flexural And Transverse Machine 150 KN



Laboratory Instrument India