



Email : sales@laboratoryinstrumentindia.com

Product Name :
Workshop Magnetic Particle Flaw Detector

Product Code :
LBNY-0004-740031



Description :

Workshop Magnetic Particle Flaw Detector

Technical Specification :

Magnetic particle flaw detector

Portable high-capacity lithium battery pack without limitation of 220V power supply. With low power consumption and small volume, charging either online or offline. It is especially proper for working in the wild and aerial area.

Adjustable illumination brightness and working power can meet different testing demand and save the setting.

New power supply can enhance the testing depth and sensitivity.

Pulse power supply can continuous work without chill time.

Combine all traditional advantages in one unit can detect surface flaws fast, clearly and accurately.

Adjustable brightness, power and range can meet different demand.

Small size with equipped bag can improve working efficiency.

Standard for welding

Working Voltage: 220 50Hz

Working Current: AC1.4A

Electrode Spacing: 0~210mm

Weight: 3.2kg

Working circle: 30 seconds on 30 seconds off

Size: 208X52X146cm

Workshop Magnetic Particle Flaw Detector, Workshop Magnetic Particle Flaw Detector Equipments, Workshop Magnetic Particle Flaw Detector Tools, Workshop Magnetic Particle Flaw Detector Tool Kits, Workshop Magnetic

Particle Flaw Detector Manufacturers, Workshop Magnetic Particle Flaw Detector Suppliers from India, China, Kenya



Laboratory Instrument India