



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Device With Central Screw

**Product Code :**  
LBNY-0005-122000126



**Description :**

Device With Central Screw can be foreseen on all models of concrete compression machines, except High Stability models. Recommended solution for machines equipped with big sized platens. Very practical to adjust the light between the compression platens of a machine, according to the height of the specimen to be tested.

**Technical Specification :**

Device With Central Screw



**Laboratory Instrument India**