



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Cube Centering Device Platen Dia

**Product Code :**  
LBNY-0005-122000128



**Description :**

The lower compression platen of the testing machine is marked with a series of concentric circles to facilitate the correct centering of the specimens. For cubes 100 and 150 mm side and cylinders dia.100 and 150 mm. Auto-Centering Device, to be used with compression machine having platen dia. 287 mm. However to grant a rapid and accurate centering of concrete cube and cylinder specimens, this Auto-Centering device is recommended.

**Technical Specification :**

Cube Centering Device Platen Dia



**Laboratory Instrument India**