



Email : sales@laboratoryinstrumentindia.com

Product Name :
Workshop UT Unit (Bend Tester)

Product Code :
LBNY-0004-740032



Description :

Workshop UT Unit (Bend Tester)

Technical Specification :

UT Unit (Bend Tester)
Standard for welding
Tabletop size for shop or off-site portability
3 plungers
Adjustable feet
Bends 1/8" plate, 1/4" plate and 3/8" plate
Height: 24.5" (622 mm)
Width: 13.75" (350 mm)
Depth: 11" (280 mm)

Workshop UT Unit (Bend Tester), Workshop UT Unit (Bend Tester) Equipments, Workshop UT Unit (Bend Tester) Tools, Workshop UT Unit (Bend Tester) Tool Kits, Workshop UT Unit (Bend Tester) Manufacturers, Workshop UT Unit (Bend Tester) Suppliers from India, China, Kenya

Laboratory Instrument India

