



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Bandsaw Butt Welding Machine

**Product Code :**  
LBNY-0005-1230004



**Description :**

**Precision Welding:**

The movable jaws run on precision ball type guides which reduce friction to almost nothing and ensure utmost welding accuracy. Minute precision control is ensured by the use of latest control devices.

**Operation:**

All the parts have been duly dimensioned; clearly arranged simple levers ensure easy operation even for unskilled workers. Sturdy in construction, Butt Welder is so designed as to facilitate reliable and easy operation even by unskilled personnel. Clamping arrangements are by quick action clamps and the alignment is automatic and perfect which is so essential in modern production methods. It is most suitable for industrial conditions and will give a longlasting dependable service.

**Technical Specification :**

Job: Carbon Steel Blades Carbon Steel Blades Bi-Metal Blades

Welding Range: 3 - 25 mm 6 - 37 mm 12-31 mm

1/8 to 1 1/4 to 1 1/2 1/2

Capacity: 3KVA (Nominal) 5KVA (Nominal) 12KVA

Input Voltage: 230V 50 Cycles 230V 50 Cycles 2.0 400V 50 Cycles

Duty Cycle: 5 5 5

---

# Laboratory Instrument India

