



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Tilting Arbour Circular Saw With Sliding Table

**Product Code :**  
LBNY-0005-12300010



**Description :**

mm/inch mm/inch  
Max circular saw diameter: 355/14 355/14  
Max.depth of cut: 114/41/2 114/4 1/2  
Arbour tilt-degree: 45 45  
Work table size: 915x785/36x31 915x785/36"x31  
Sliding table size: 710x400/28/15 3/4 -  
Sliding table movement: 1000/39 -  
Size of dust collection port-dia: 100/4 100/4  
Circular saw speed (rpm): 4000 4000  
Power required hp/rpm: 3/1440 3/1440

**Technical Specification :**

Tilting Arbour Circular Saw With Sliding Table



**Laboratory Instrument India**

