



Email : sales@laboratoryinstrumentindia.com

Product Name :
Chain and Chisel Mortiser

Product Code :
LBNY-0005-12300012



Description :

mm/inch
Table size: 455x 150/18x6
Max.height of workpiece on table: 254/10
Max.height of door or frame: 1000/39
Max.depth of cut-chain mortiser: 140/51/2
Max.depth of cut-chisel mortiser: 90/31/2
Longitudinal movement of head/table: 305/12
Table tilt-degree: 45
Size of dust collection port-dia: 100/4
Air pressure required (kg/cm²): 5 to 7
Electric motor (hp/rpm): 2/1440

Technical Specification :

Chain and Chisel Mortiser



Laboratory Instrument India

