



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Chain and Chisel Mortiser

**Product Code :**  
LBNY-0005-12300012



**Description :**

mm/inch  
Table size: 455x 150/18x6  
Max.height of workpiece on table: 254/10  
Max.height of door or frame: 1000/39  
Max.depth of cut-chain mortiser: 140/51/2  
Max.depth of cut-chisel mortiser: 90/31/2  
Longitudinal movement of head/table: 305/12  
Table tilt-degree: 45  
Size of dust collection port-dia: 100/4  
Air pressure required (kg/cm<sup>2</sup>): 5 to 7  
Electric motor (hp/rpm): 2/1440

**Technical Specification :**

Chain and Chisel Mortiser



**Laboratory Instrument India**

