



Email : sales@laboratoryinstrumentindia.com

Product Name :
Wood Saw

Product Code :
LBNY-0005-12300013



Description :

mm/inch mm/inch mm/inch mm/inch
Diameter of saw wheels: 305/12 455/18 305/12 455/18
Throat depth: 280/11 430/17 280/11 430/17
Max. cutting capacity: 150/6 254/10 150/6 254/10
Work table size: 380x380/15x15 535x585/21x23 355x380/14x15 508x560/20x22
Table tilt degree: 45 45 45 45
Width of saw wheels: 38/1 1/2 50/2 38/1 1/2 50/2
Speed of saw wheels (rpm): 500/750 500/750 500/750 500/750
Power required hp/rpm: 1.5/1440 2/1440 1.5/1440 2/1440

Technical Specification :

Wood Saw



Laboratory Instrument India