




Email : sales@laboratoryinstrumentindia.com

Product Name :
Single Column Vertical Lathe

Product Code :
LBNY-0005-12300019



Description :

Max. Cutting Dia. For Vertical Tool Post: 1000mm 1250mm 1600mm 2000 mm
Max. Cutting Dia. For sideway tool past: 1000mm 1100 mm 1400mm 1800mm
Table dia.: 900mm 1000mm 1400mm 1800mm
Max. Work height: 800mm 1000mm 1000mm 1250mm
Max. work load: 2000 Kg. 3200Kg. 5000Kg. 8000Kg.
Table speed range: 16/8-250rpm 16/6.3-200rpm 16/5-160rpm 16/4-125rpm
Feed range for both tool past: 0.9-103 0.8-86 0.8-86 0.8-86
Horizontal travel for vertical tool past: 670mm 700mm 900mm 1110mm
Vertical travel for vertical tool past: 650mm 650mm 800mm 800mm
Horizontal travel for sideway tool past: 500mm 500mm 650mm 630mm
Vertical travel for vertical tool past: 810mm 900mm 900mm 1110mm
Ram travel: 580mm 650mm 650mm 650mm
Main motor power: 22kW 22kW 22kW 30kW
Machine Weight: 7100Kg. 8500Kg. 13000 Kg. 17500Kg.

Technical Specification :

Single Column Vertical Lathe

Laboratory Instrument India

