



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Hand Operated Bead Bending

**Product Code :**  
LBNY-0005-12300025



**Description :**

**Features:**

Special steel adjustable.  
4-sets standard rollers.  
Heavy and solid cast iron construction.

**Technical Specification :**

Thickness max: 1.2mm  
Cylinder length: 140mm  
Throat Depth: 200mm  
Net weight: 48kg  
Packing size: 670x320x600mm



**Laboratory Instrument India**