



Email : sales@laboratoryinstrumentindia.com

Product Name :
Double Column Vertical Lathe

Product Code :
LBNY-0005-12300026



Description :

Max. Cutting Dia. For Vertical Tool Post: 2500mm 3150mm 4000mm 5000mm
Working table diameter: 2250mm 2830mm 3200mm 4000/4500mm
Max. height of work-piece: 1600/2000mm 1600/2000mm 2000/2500mm 2500/3150mm
Max. weight of work-piece: 10 to 20Ton 10 to 20Ton 20Ton 32/50Ton
Working table range of rotation speed: 2 to 63rpm 2 to 63rpm 2 to 63rpm 0.6-25.4rpm
Working table step of rotation speed: 16 steps 16 steps 16 steps stepless
Horizontal travel of tool post: 1400mm 1650mm 2150mm 2550mm
Crossbeam travel: 1250/1650mm 1250/1650mm 1650/2150mm 2150mm
Power of main motor: 55kw 55kw 55kw 75kw
Weight of machine(approx.): 41/43 Ton 43/45 Ton 45/46 Ton 57/80 Ton

Technical Specification :

Double Column Vertical Lathe



Laboratory Instrument India

