



Email : sales@laboratoryinstrumentindia.com

Product Name :
Power Operated Bead Bending

Product Code :
LBNY-0005-12300030



Description :

Thickness Max.: 1.2mm 2.5mm 4.0mm
Cylinder Length: 140mm 120mm 165mm
Throat Depth: 200mm 280mm 300mm
Cylinder Speed: 32rpm 20rpm 15rpm
Motor Power: 1 HP 2 HP 3 HP
Net Weight: 120kg 235kg 341kg
Packing Size: 1100x480x1480mm 1470x580x1350mm 1650x580x1400mm

Technical Specification :

Power Operated Bead Bending



Laboratory Instrument India