



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hydraulic Hand and Motorised Bending Machine

Product Code :
LBNY-0005-12300031



Description :

Capacity: 3/8, 1 1/2 to 2 1/2 to 3
Pressure: 10 Ton 15 Ton 15 Ton
Former Size: 3/8, 1/2, 1/2, 3/4, 1/2, 3/4
3/4, & 1 1/4, 1.1/4, T, 1.1/4,
1.1/2, & 2 1/2, & 2
2.1/2, & 3

Technical Specification :

Hydraulic Hand and Motorised Bending Machine



Laboratory Instrument India