



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Hydraulic Iron Worker

**Product Code :**  
LBNY-0005-12300033



**Description :**

**Silent Features:**

Unique compact design by using ultrasonically tested steelplates.  
International brands hydraulic & electrical components ensure trouble free operation.  
Hydraulic Iron Worker is designed to give an effective solution to all the cutting and punching problems of flat and Section with satisfactory work quality from robust and compact machine.  
Direct shearing avoids any waste of materials.  
Five work station catering 7 different operation. Parallel operation of punching station with other work station gives higher productivity.

**Technical Specification :**

Hydraulic Iron Worker



**Laboratory Instrument India**