



Email : sales@laboratoryinstrumentindia.com

Product Name :
Workshop Computer Monitor (Flat Touch Screen)

Product Code :
LBNY-0004-750003



Description :

Workshop Computer Monitor (Flat Touch Screen)

Technical Specification :

Computer monitor (Flat touch screen)

Built-in screen that reduces potentially harmful blue light emissions while delivering excellent color accuracy. Color consistency comes standard: See consistent, vibrant colors across a wide viewing angle enabled by In-Plane Switching (IPS) technology. With 99% sRGB coverage, you get precise color right out of the box.

Sleek, stylish and organized: Hide cords in the monitor riser to keep things organized. This sleek monitor has a three-sided ultrathin bezel, small base and an improved, easy-to-use cable management system.

Comfort is key: Tilt, swivel, pivot and adjust the height of your monitor (150mm height adjustment range) for extra comfort during work. Designed with you in mind: Snap on the Slim Soundbar and elevate your listening experience.

Flat touch Screen

size 22inches

Ram 6gb

Hard disk 1TB

Workshop Computer Monitor (Flat Touch Screen), Workshop Computer Monitor (Flat Touch Screen) Equipments, Workshop Computer Monitor (Flat Touch Screen) Tools, Workshop Computer Monitor (Flat Touch Screen) Tool Kits, Workshop Computer Monitor (Flat Touch Screen) Manufacturers, Workshop Computer Monitor (Flat Touch Screen) Suppliers from India, China, Kenya



Laboratory Instrument India