



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Hand Operated Folding Machine

**Product Code :**  
LBNY-0005-12300035



**Description :**

**Features:**

Compact High Quality Engineering.  
Ensure Minimum Maintenance, Together With Simple and Safe Operation.  
A Fully Welded Structure.  
With Foot Control.  
It Is Easy For Operation And Relax Hands.

**Technical Specification :**

Maximum Working Length: 1020mm 1270mm 1520mm 2020mm 2520mm  
Maximum Sheet Thickness: 2.5mm 2.0mm 1.5mm 1.2mm 1.0mm  
Maximum Clamping Bar Lift: 47mm 47mm 47mm 47mm 47mm  
Folding Angle: 0 - 135 0 - 135 0 - 135 0 - 135 0 - 135  
Net Weight: 285Kgs. 320Kgs. 385Kgs. 490Kgs. 590Kgs.



**Laboratory Instrument India**

