



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hand Operated Folding Machine

Product Code :
LBNY-0005-12300035



Description :

Features:

Compact High Quality Engineering.
Ensure Minimum Maintenance, Together With Simple and Safe Operation.
A Fully Welded Structure.
With Foot Control.
It Is Easy For Operation And Relax Hands.

Technical Specification :

Maximum Working Length: 1020mm 1270mm 1520mm 2020mm 2520mm
Maximum Sheet Thickness: 2.5mm 2.0mm 1.5mm 1.2mm 1.0mm
Maximum Clamping Bar Lift: 47mm 47mm 47mm 47mm 47mm
Folding Angle: 0 - 135 0 - 135 0 - 135 0 - 135
Net Weight: 285Kgs. 320Kgs. 385Kgs. 490Kgs. 590Kgs.



Laboratory Instrument India

