



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
H Type Power Operated Hydraulic Press

**Product Code :**  
LBNY-0005-12300039



**Description :**

**Features:**

**Structure:**

The Presses have good rigid structure for accurate works.

**Power Pack Unit:**

Our Presses are powered by specially designed Power Pack Unit having Oil Breather, Oil Level Indicator, Electric Motor, Hand Lever operated Direction Control Valve, Suction Strainer, Gear pump, Flexible Coupling, Pressure Gauge, Oil Tank (without oil), Return Line Filter, Internal & External Pipe Line etc.

**Hydraulic Cylinder:**

All chevron seals & bronze strip used for leak proof performance. Ram will be made from material & surface ground hardened for long life of seals & corrosion resistance. Our Presses are powered by specially designed Double Acting Hydraulic Cylinders. Cylinders manufactured out of quality steel. Cylinder bores will be honed to mirror finish or long life of piston seals.

**Technical Specification :**

H Type Power Operated Hydraulic Press



**Laboratory Instrument India**

