



Email : sales@laboratoryinstrumentindia.com

Product Name :
H Type Hand Operated Hydraulic Press

Product Code :
LBNY-0005-12300041



Description :

Fast and safe adjustment of the press table height using a braked winch. Hand operated workshop press. Hard chrome piston and very long stroke of 300mm. Available from 40 to 100 tons. Foot operated air approach feature.

Machine Features:

Fast and safe adjustment of the press table height using a braked winch
Movable cylinder
Foot operated air approach feature
Fast air return
Wide insertion width
Hard chrome piston and very long stroke of 300mm

Technical Specification :

Maximum daylight: 750mm
Width between supports: 840mm
Maximum (ton): 50 tonne
mm/stroke: 6,0 / 2,0
Table width: 200mm
Ram diameter: 70mm
Cylinder ram stroke: 300mm
Cylinder diameter: 125mm
Height overall: 1990mm
Width overall: 1350mm

Depth overall: 630mm
Weight: 370kg.



Laboratory Instrument India