



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
H Type Hand Operated Hydraulic Press

**Product Code :**  
LBNY-0005-12300041



**Description :**

Fast and safe adjustment of the press table height using a braked winch. Hand operated workshop press. Hard chrome piston and very long stroke of 300mm. Available from 40 to 100 tons. Foot operated air approach feature.

**Machine Features:**

Fast and safe adjustment of the press table height using a braked winch  
Movable cylinder  
Foot operated air approach feature  
Fast air return  
Wide insertion width  
Hard chrome piston and very long stroke of 300mm

**Technical Specification :**

Maximum daylight: 750mm  
Width between supports: 840mm  
Maximum (ton): 50 tonne  
mm/stroke: 6,0 / 2,0  
Table width: 200mm  
Ram diameter: 70mm  
Cylinder ram stroke: 300mm  
Cylinder diameter: 125mm  
Height overall: 1990mm  
Width overall: 1350mm

---

Depth overall: 630mm  
Weight: 370kg.



**Laboratory Instrument India**