



Email : sales@laboratoryinstrumentindia.com

Product Name :
C Type Power Operated Hydraulic Press

Product Code :
LBNY-0005-12300042



Description :

Structure:

The Presses have good rigid structure for accurate works.

Hydraulic Cylinder:

All chevronseals & bronze strip used for leak proof performance. Ram will be made from material & surface ground hardened for long life of seals & corrosion resistance. Our Presses are powered by specially designed Double Acting Hydraulic Cylinders. Cylinders manufactured out of quality steel. Cylinder Bores will be honed to mirror finish or long life of piston seals.

Power Pack Unit:

Our Presses are powered by specially designed Power Pack Unit having Oil Breather, Gear Pump, Flexible Coupling, Electric Moto, Hand Lever operated Direction Control Valve, Oil Level Indicator, Suction Straine, Pressure Gauge, Oil Tank, Return Line Filter, Internal & External Pipe Line etc.

Technical Specification :

C Type Power Operated Hydraulic Press

Laboratory Instrument India

