



Email : sales@laboratoryinstrumentindia.com

Product Name :
Band Saw Blade Sharpener

Product Code :
LBNY-0005-12300045



Description :

It is suitable for grinding the bandsaw that the teeth are same. Automatically dial tooth, easy to operate high precision. Make the old band saw reuse and save cost. It can grind various thin wood saw blade, double metal sawblade, carbide steel saw blade and so on. It is a new multifunctional patent machine that can grind changeable teeth (tip) of band saw blade, such as 2/3, 1 4/6 etc.

Technical Specification :

Saw blade width: 16-45mm
Teeth space: 3-30mm
Feed speed: 40-60teeth/min
Grinding wheel: 150mm
Power/Speed: 380V/220V/550W/2800RPM
Dimension/Weight: 60x50x 100cm/80kg



Laboratory Instrument India