



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Band Saw Blade Sharpener

**Product Code :**  
LBNY-0005-12300045



**Description :**

It is suitable for grinding the bandsaw that the teeth are same. Automatically dial tooth, easy to operate high precision. Make the old band saw reuse and save cost. It can grind various thin wood saw blade, double metal sawblade, carbide steel saw blade and so on. It is a new multifunctional patent machine that can grind changeable teeth (tip) of band saw blade, such as 2/3, 1 4/6 etc.

**Technical Specification :**

Saw blade width: 16-45mm  
Teeth space: 3-30mm  
Feed speed: 40-60teeth/min  
Grinding wheel: 150mm  
Power/Speed: 380V/220V/550W/2800RPM  
Dimension/Weight: 60x50x 100cm/80kg



**Laboratory Instrument India**