



Email : sales@laboratoryinstrumentindia.com

Product Name :
Saw Blade Sharpener

Product Code :
LBNY-0005-12300046



Description :

It can grind the teeth. And it also can keep balance and efficiently. Without model plate, can grind the teeth automatically. Adjustable speed, it is fast and accurate.

Technical Specification :

Saw blade Dia.: 350-500 60-350
Speed: 45teeth/min
Voltage: 220V/380V 220V
Power: 550W/60W 250W/60W
Dimension: 60x60x 100cm 62x62x72cm
Weight: 75kg 56kg



Laboratory Instrument India