



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Arc Angle Chamfer

**Product Code :**  
LBNY-0005-12300047



**Description :**

Applying tonon-heat treatment workpiece for such as various mould, part of machines chamfer. Standard square alloy bit , easy to replace. Or you can use old bits that used inmilling machine The dump railadopt good heat-treatment steel plate, which will abide rough use This is originative to chamfer the R angle of workpiece.

**Technical Specification :**

Chamfering height: R1-R4  
Power: 380V 250W  
Speed: 2800rpm  
Dimension: 53x51 x24cm  
Weight: 28.8kg/32.8k



**Laboratory Instrument India**