



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
High Speed Sliding Chamfer

**Product Code :**  
LBNY-0005-12300048



**Description :**

The machine adopt Swedish bearing and imported digital cuttings.  
Easy to adjust, use linear motion guide, can protect the surface of the products.  
Also can fix the table, reduce the cutting resistance and secondary burrs.

**Technical Specification :**

Chamfering height: 0-1.5mm  
Chamfering angle: 45°  
Power: 380V 250W  
Speed: 8500rpm  
Dimension: 60 x 29 x 37cm  
Weight: 40kg



**Laboratory Instrument India**