



Email : sales@laboratoryinstrumentindia.com

Product Name :
Portable Chamfer

Product Code :
LBNY-0005-12300050



www.laboratoryinstrumentindia.com

Description :

The dump rail adopt good heat-treatment steel plate, which will abide rough use. Suitable for any length of non-heat treatment work piece and any other parts which are hard to move. It is easy to operate, take less time. Standard square alloy bit, easy to replace. Or you can use old bits that used in milling machine.

Technical Specification :

Chamfering height: 0-9mm
Chamfering angle: 15 -45 deg.
Power: 380V 250W
Speed: 2800rpm
Dimension: 38 x 24 x 24,5cm
Weight: 12kg



Laboratory Instrument India